


INSTRUCTIONS FOR USE

ToFFEE Bluetooth Switch

Pack Contents

Carefully check the contents of the box, which are:

- ToFFEE Bluetooth Switch unit
- Decal set
- USB Charging Cable
- These instructions for use



Product Description

ToFFEE BT Switch is a unique and highly versatile proximity switch designed for disabled users. It works with most Bluetooth enabled devices, including PCs, Macs, Chromebooks, iPad/iPhone and Android phones and tablets. Benefitting from ten user-selected detection ranges and output modes for a wide range of use-cases, **ToFFEE BT** is ideal not only for typical switch applications but also for situations where conventional switches fall short.

Unlike other proximity switches that require contact or very close proximity, **ToFFEE BT** uses advanced Time-of-Flight (ToF) technology to detect any opaque object, regardless of colour, material, or surface texture. This means it works reliably with any skin tone, and is unaffected by clothing, hair, or coverings.

Optional audible and visual feedback means it can be made to sound just like a conventional switch. The switch is powered by an internal rechargeable lithium battery which gives around 25 hours of use between charges. Programmable power settings allow the user to select if and when the unit goes to sleep due to inactivity.

ToFFEE BT may be mounted using industry-standard brackets, allowing it to be operated by almost any part of the body such as the user's hand, arm, knee, head or even by puffing out a cheek.

Like all Pretorian products, we've packed in as much functionality as we can, so please take the time to read these instructions to get the most from your **ToFFEE BT** switch. Your attention is drawn to the Safety Instructions in particular, which are overleaf.


IMPORTANT SAFETY INSTRUCTIONS



- **Read and pay heed to these instructions for use and retain for future reference.**
- **Ensure you are familiar with the devices you use ToFFEE BT Switch with by reading their instructions for use and/or using them separately to begin with.**
- **This product contains a laser light source, which is designed to meet Class 1 laser safety limits under all reasonably foreseeable conditions including single faults, in compliance with IEC 60825-1:2014. This product is not considered 'child appealing' within the context of Clause 5 of IEC60825-1:2014. Nevertheless, you should ensure that users do not stare into the beam at very close range. If they are unable to follow such instruction, or the unit cannot be positioned in such a way as to make it impossible, this constitutes a contra-indication for use.**
- **Whilst ToFFEE BT Switch uses state of the art optical technology, it is still possible for it to be activated unintentionally due to external light sources or electrical interference. Do not use it in applications where unintentional operation may lead to safety concerns.**
- **The cover glass [1] should be kept free of debris and fingerprints for reliable operation. Refer to the cleaning guide in these instructions for use.**
- **There are no user serviceable parts in the unit. Do not disassemble. Refer to an authorised distributor or the manufacturer for all repairs.**
- **Do not use this unit near water.**
- **This unit is intended for indoor use only.**
- **This unit is not for use as part of a safety critical, mission critical or life-sustaining system.**
- **ToFFEE BT Switch contains a lithium-ion battery. Never dispose of a battery in fire.**

Features

- **Choice of ten ranges: 0 - 25mm up to 0 - 500mm**
- **Choice of 26 output functions, including space, enter, media, camera etc.**
- **Operates with any skin tone or opaque material**
- **Multiple output options for a wide range of use-cases**
- **Ignore Repeats setting**
- **Choice of audible and visual feedback, both, or none.**
- **Rechargeable lithium battery with supplied charger cable**
- **Selectable power-down options**
- **Threaded mounting holes.**



Preparing your ToFFEE BT switch for use

First of all, choose the colour you would like **ToFFEE BT** to be and attach the self-adhesive decal to the top surface, positioning it carefully to ensure good alignment with the recess in the mouldings.

Make sure your **ToFFEE BT** is fully charged using the cable provided. Simply insert the USB-C plug into the corresponding socket [4] and plug the other end into a computer or charger socket. Once charging begins, the charging LED [8] lights green. When this LED goes out, charging is complete.

Note that if you should lose your charger cable, replacements can be purchased either from a Pretorian reseller or any electronics or on-line retailer. You should request a 'USB-A to USB-C cable', (or USB-C to USB-C if your charger also uses USB-C).

Once charged, you may connect it and make the required settings so that it behaves as you need it to. Please refer to the following sections for instructions on how to configure each of the settings.

Connecting to your device

When new, or when **ToFFEE BT** has been put to sleep manually, wake it up by pressing the Function/Modify/Power button [7]*.

Go to the Bluetooth menu on your device (on iDevices, Settings → Bluetooth) and first of all ensure Bluetooth is turned on (on iDevices, using the slider at the top of the screen).

After a few seconds the unit should appear as a 'discoverable' device. It will appear as something similar to:

ToFFEE BT-V164.0-ABC1

Tap or click on the name and the pairing process will begin. Typically, it takes around 5 seconds to connect, after which the device will state that it is 'Connected'.

***ToFFEE BT** is shipped in Sleep mode- please refer to the Power Settings and Manual Sleep sections for more details.

Range Setting

Range is the farthest distance from the cover glass that the unit can detect an object. In all cases, the unit will detect objects from zero (i.e. in contact with the unit) up to and including the range setting.

To see the current Range, press the Range/Learn button [6] briefly to bring on the display. The corresponding Range settings are given in Table 1. If you wish to change the setting, press it repeatedly (but briefly) until you see the setting you need, and then wait for the display to be extinguished. You can gauge the effect of changes to the Range setting immediately, since both audible and visual feedback remain active (when enabled). However, the unit does not output anything over Bluetooth while setting the Range. Note that Range settings are approximate.

Range Setting	Range (millimetres)	Range (inches)
0	0-25 mm	0-1"
1	0-50 mm	0-2"
2	0-75 mm	0-3"
3	0-100 mm	0-4"
4	0-150 mm (default)	0-6" (default)
5	0-200 mm	0-8"
6	0-250 mm	0-10"
7	0-300 mm	0-12"
8	0-400 mm	0-16"
9	0-500 mm	0-20"

Table 1:
Range settings

Function Setting

Function sets the character or command that is sent to the computer whenever the switch is activated. There are 26 settings, which are a carefully curated list of settings that are appropriate to the hardware and software that is commonly used by disabled users.

To see the current Function, press the Function/Modify button [7] briefly to bring on the display. The corresponding Function settings are given in Table 2. If you wish to change the setting, press the button repeatedly (but briefly) until you see the setting you need, and then wait for the display to be extinguished.

No Bluetooth outputs occur while the Function is being set but become 'live' again as soon as the display is extinguished. They are also saved in the unit's memory for future use.

Function Setting	Function	Type	Allowable Output Modes
0	Numeral 1	Keyboard	All
1	Numeral 2	Keyboard	All
2	Numeral 3	Keyboard	All
3	F7	Keyboard (e.g. for Clicker)	All
4	F8	Keyboard (e.g. for Clicker)	All
5	Space (Default)	Keyboard (widely used for switch adapted app access)	All
6	Enter	Keyboard (widely used for switch adapted app access)	All
7	~ 1	Keyboard	Forces One-Shot
8	~ 3	Keyboard	Forces One-Shot
9	Up Arrow	Keyboard	All
A	Down Arrow	Keyboard	All
B	Left Arrow	Keyboard	All
C	Right Arrow	Keyboard	All
D	Keyboard Up/Down	Device	Forces One-Shot
E	Play/Pause	Media (for music and video)	Forces One-Shot
F	Track Forward/ Seek Forward	Media	Double-click not permitted- forces Direct
G	Track Back/ Seek Backwards	Media	Double-click not permitted- forces Direct
H	Volume Up	Media (also takes photos/video)	All
J	Volume Down	Media (also takes photos/video)	All
L	Mute	Media	Forces One-Shot
N	Timed Play 10s	Media (often used to encourage reluctant switch users)	Forces One-Shot
P	Timed Play 30s	Media (often used to encourage reluctant switch users)	Forces One-Shot
R	Home	Device (equivalent to swipe up from bottom on iPad)	All
T	Enter/Home	Combo (Enter on short press, Home on long press)	Forces Direct
U	Left Click	Mouse (can be combined with Latch, Double Click or One-Shot)	All
Y	Right Click	Mouse (can be combined with Latch, Double Click or One-Shot)	All

Table 2: Function settings

Users who are familiar with APPLICATOR and iSwitch/ iSwitch50 from Pretorian will be immediately familiar with most of this list.

Play/Pause, Track Forward and Track Back may be used with music and media apps such as Apple Music, Spotify and YouTube. Note that keeping Track Forward and Track Back activated enables the user to seek within a track. If that is undesirable functionality (for example because the user takes time to remove their hand from the beam), enable One-Shot mode - see below. To make sure that media commands are directed to the correct app, it's best to make sure that only one such app is open at once.

Volume Up and Volume Down may also be used to take photos and start/stop video recording in the camera app, which can be a great way of engaging first-time users.

Timed Play 10s and Timed Play 30s allow you to create a 'reward', for example when completing a task. It is particularly useful in helping to engage reluctant switch users. Because this setting uses the 'Play/Pause' command, which is a 'toggle action', it is important that the device is paused (not playing) before you press the switch, otherwise it will pause for a timed period instead of playing.

Left Click and Right Click may frequently need to be used in combination with the Output Mode settings to give Drag, Double-Click or One-Shot (momentary click) functions – see below.

ToFFEE BT will often be used with iOS Switch Control. When using it in this way, choose Functions that give single, 'normal' keyboard outputs such as Numerals 1-3, Space and Enter. Avoid using -1 and -3, since they send multiple keystrokes. Additionally, media functions such as Play/Pause are not recognised as valid keystrokes by Switch Control.

The Enter/Home setting can be especially useful for iOS Switch Control users who make use of the auto-scanning feature, since it allows everything to be controlled with a single switch and the app then closed afterwards.

Alternatively, some users may use **ToFFEE BT** with iOS Assistive Touch, in which case Left Click is the correct Function. Assistive Touch does not currently recognise Right Click.

The Range and Function settings may often need to be changed, so they have been made easily accessible from the keypad. The remaining settings usually need to be changed less often and so are bundled together into the Learn Mode:

Learn Mode

In Learn Mode, you can set the following

- Output Mode
- Ignore Repeats time
- Feedback type
- Power settings

To enter Learn Mode, press and hold the Range/Learn button [6] until you hear a warbling beep.

Output Mode

The Output mode is the first setting to appear on the display, the settings being as shown in Table 3:

Which setting is appropriate very much depends on the use-case. For example, Double-Click is particularly useful in combination with Left Click function for computer access, since it can emulate a double mouse click for those users who lack to dexterity to click twice in short succession. Latching is equally useful as it can give a Drag feature.

One-shot is useful with a variety of Functions for users whose movements are laboured and who will tend to linger over the switch for longer than the period for which they wish to operate it. Additionally, it can be used to good effect in preventing multiple operation of **ToFFEE BT** when an object is introduced to the very edge of the sensor zone and might otherwise cause the switch to 'chatter'.

Delayed Drag is a useful blend of Direct and Latched. If an object passes through the sensor zone and back out in less than 1.5 seconds, it operates in Direct mode. If the object remains in the beam for more than 1.5 seconds, the switch latches on. If audible feedback in on (see later), the unit will click a second time to indicate that the latch is now on. The output can be de-latched by introducing an object to the sensor zone again, even for a short period.

Whilst having separate Function and Output Mode settings give maximum flexibility, there are some combinations that don't make sense, or may lead to operator confusion. To prevent this, the unit sometimes overrides the Output Mode setting. For example, the combination of Play/Pause and Double-Click would lead to the unit playing and then almost immediately stopping again. To prevent this, selecting Play/Pause automatically forces the unit to work in One-Shot without needing to change the Output Mode. If you navigate away from Play/Pause to a function which can accept the actual Output Mode setting, it will automatically change back. This is summarised in the last column of Figure 2.

A	Direct/momentary (default)
B	Latching
C	One-Shot
D	Double click
E	Delayed Drag

Table 3:
Output Mode settings

To make changes in this menu, press the Function/Modify button [7] repeatedly until you see the setting you require on the display. Then, either let the display extinguish or press the Range/Learn button [6] to go to the next Learn Mode menu.

Ignore Repeats time

This is the time for which the switch will ignore objects after it was last released. It may be used to prevent accidental secondary operation of the switch, for example because the user has a tremor. By default, it is set to zero, i.e. the unit will operate again immediately after the switch is released.

When Ignore Repeats is set to a non-zero value, the red indicator LEDs [3] give a visual indication that objects will be ignored following the last switch operation. Once these LEDs are extinguished, the switch will respond to objects once again. The settings are as follows:

F	None (default)
G	0.5 sec
H	1.0 sec
J	1.5 sec
L	2.0 sec

Table 4:
Ignore Repeats
settings

To make changes in this menu, press the Function/Modify button [7] repeatedly until you see the setting you require on the display. Then, either let the display extinguish or press the Range/Learn button [6] again to go to the next Learn Mode menu:

Feedback type

The unit is capable of both visual and audible feedback. Visual feedback is achieved using four green LEDs [2] on the top of the switch. Audible feedback has been designed to sound just like a normal, mechanical switch. Each can be enabled and disabled independently using this menu. The settings are given in Table 5.

N	None
P	Audible only
R	Visual only
S	Audible and Visual (default)

Table 5:
Feedback settings

To make changes to this setting, press the Function/Modify button [7] repeatedly until you see the setting you require on the display. Then, either let the display extinguish or press the Range/Learn button [6] again to go to the next Learn Mode menu:

Power Settings

Three power regimes are available. Which to use depends on how tolerant the user is of the inconvenience of the unit going to sleep. Setting T is perhaps the best compromise of battery life and convenience, but it does require that the Power button be pressed after very long periods of inactivity, and not all users will be able to do so. That situation is likely exacerbated when the unit is mounted, as this will make access to the Power button more difficult.

When set to Power Setting T, the unit goes into Snooze mode after 30 minutes of inactivity. It emits two short beeps to indicate it is doing so, regardless of whether audible feedback is turned on. If there is no further activity for another 2.5 hours, it emits four beeps and then goes to Sleep.

During Snooze, the unit is disconnected from Bluetooth. It will wake up when a suitable object is introduced at a distance between zero and the current Range setting and the Bluetooth connection will automatically be restored. Power savings are also achieved by 'looking' for an object less often, so you will need to keep the object in range for longer than usual to make sure it is recognised. As it awakes, the unit beeps twice, regardless of the audible feedback setting.

During Sleep, the unit does not respond to objects at all. The only way to wake the unit is to press the Function/Modify/Power button [7]. Battery life in Sleep mode is many months on a single charge. As with Snooze mode, Sleep disconnects the unit from Bluetooth.

Power Setting U allows Snooze to be entered after 30 minutes (just like setting T) but the unit will never go to Sleep. This is a good choice for users who cannot press the Power button due to their disability, or because the unit is mounted. To wake the unit from Snooze, introduce an object at a distance between zero and the Range setting. As with setting T, you may need to keep it in range for slightly longer when waking from Snooze. The unit beeps twice as it awakes. Battery life in Snooze mode is about a week.

Power Setting Y keeps the unit awake indefinitely. It is recommended only when the unit is permanently powered via the USB socket. Battery life when Awake is about 25 hours, or indefinite if plugged into a supply.

T	Snooze after 30 minutes, Sleep after 3 hours (default)
U	Snooze after 30 minutes. Never Sleeps.
Y	Permanently Awake (recommended only when powered via cable)

Table 6: Power settings

To make changes to this setting, press the Function/Modify button [7] repeatedly until you see the setting you require on the display. Then, either let the display extinguish or press the Range/Learn button [6] again. This is the last of the Learn Mode menus, so a further press will circle back to the Output Mode setting.

On-Screen Keyboard on iDevices

Because **ToFFEE BT** appears to an iDevice as a keyboard, it automatically turns off the on-screen keyboard. This can cause inconvenience when using apps that require typed input.

To overcome this, the unit allows you to deploy the on-screen keyboard manually at any time. Beginning with the display [2] off, simply press and hold the Function/Modify button [7]. After 1.5 seconds, the bottom segment of the display will light. Now release the button and the on-screen keyboard will be deployed automatically. To turn it off again, press and hold the Function/Modify button once again until the bottom segment lights, then release.

Mode setting 'D' allows **ToFFEE BT** to be configured to deploy the on-screen keyboard using the main switch [1], although it is usually only useful when multiple Bluetooth switches are present and paired with the device.

iDevices remember your on-screen keyboard preference so there's no need to press to deploy every time.

Note that iDevices only allows the on-screen keyboard to be deployed when a text input box is selected.

Manual Sleep

As well as the power settings described above, there is also the option to put the unit to sleep manually. This can be particularly useful when the unit is to be transported (for example in a moving vehicle), to prevent unintentional activations of the switch waking the unit and reducing battery charge.

Simply press and hold the Function/Modify/Power button [7] until all the segments on the display light (about three seconds), then release*. **ToFFEE BT** goes into Sleep Mode and will not wake until the Function/Modify/Power button [7] is pressed once again.

* Note that after 1.5 seconds of pressing and holding the Mode button, the bottom segment of the display [2] lights to show that it has been pressed for long enough to send the Keyboard Up/Down command. If you intend to place the unit into Sleep Mode, do not release the button at this stage. Instead, continue pressing until ALL the display segments light, then release.

Battery Life and Battery Charging

A fully charged battery gives approximately 25 hours of operation, depending on the amount of use. When the battery is getting low, the red LEDs [3] begin to blink red. This is an indication that you should charge the battery soon. Plug the charging cable into the charging socket [4] and then into a charger or computer*.

While charging, the charging LED [8] will be illuminated green. Once charging is complete (a matter of a around 2-3 hours if fully discharged) the charging LED will be extinguished. You can then unplug the cable.

Note that you can continue using **ToFFEE BT** while it is charging.

* **ToFFEE BT** plugs into a computer USB port for charging purposes only - it does not give a functional connection this way.

Mounting

ToFFEE BT may be mounted where necessary using the two threaded mounting holes on the underside of the unit. The distance between them is 31.75mm (1.25"), which is frequently referred to by mounting manufacturers as the 'Buddy button' system. Custom mounting plates are available for ToFFEE Wired Switch and ToFFEE BT from Rehadapt and their distributors. Visit <https://rehadapt.com/>

Please note that the threaded holes require 3mm metric screws. Maximum length inside the unit is 8mm. In North America, please refer to a Pretorian reseller for availability of suitable metric screws.

The cover glass [1] must be kept clean to ensure correct functionality. Cleaning may be achieved with a slightly dampened soft cloth. Do not use any detergents, cleaning chemicals or solvents, the use of which will void the warranty.

Notes on using ToFFEE BT

ToFFEE BT uses invisible laser light to detect the presence of an object. The small black ‘cover glass’ (actually plastic) [1] on the top of the switch is where the light is emitted from and received back again. The electronics inside the unit measures the time it takes for the light to travel to the object and be reflected back again (its Time of Flight – hence ToFFEE).

Although not visible light, it is a laser, and it is important not allow anyone to stare directly into the beam at close range. If the user cannot be so instructed, or prevented from doing so, this constitutes a contra-indication for its use.

Although visible lasers are usually associated with straight beams, the laser in ToFFEE BT does not operate in that fashion. Instead, it propagates as a cone of light with an 18° ‘field of view’ from the centre of the cover glass, as shown in Figure 1. This means that as the object moves further away from the unit, the breadth of the sensor zone increases. Table 7 gives the breadth of the sensor zone at the maximum distance for each Range setting:




Figure 1 – Sensor zone

Range Setting	Maximum distance (L in Figure 1)	Diameter of sensor zone (D in Figure 1)
0	25 mm	8 mm
1	50 mm	16 mm
2	75 mm	24 mm
3	100 mm	32 mm
4	150 mm	48 mm
5	200 mm	64 mm
6	250 mm	80 mm
7	300 mm	96 mm
8	400 mm	128 mm
9	500 mm	160 mm

Table 7:
Sensor zone
maximum
diameter

You can consider the sensor zone as a cone emanating from the cover glass with a flat top at the Range distance and maximum diameter given in Table 7. An object anywhere in the sensor zone will activate the switch. You can therefore use three generalised methods to operate **ToFFEE BT**:

- A 'pass-through' where the object (perhaps a hand) is passed horizontally through the beam and out the other side. When using this technique, it is best to make sure the Range is around twice the nominal distance of the object from the cover glass. See Figure 2a.
- A 'wave' in which the object enters the beam and then back out the same side. Again, it's best to make the range approximately twice the nominal distance of the object. See Figure 2b.
- A 'top-down' activation where the object enters the sensor zone from above the flat top of the cone and leaves by backing away from it. Note that as the object enters the sensor zone, the switch activates when the object reaches the flat top of the cone. As the object is removed, the switch de-activates at a greater distance than this to prevent switch chatter. This is known as hysteresis and is shown diagrammatically in Figure 2c and 2d.

Which of these techniques is best very much depends on the user's abilities, but is also affected by the way in which **ToFFEE BT** is positioned. Of course, the techniques described above can be used interchangeably.

Users who lack fine motor skills or are subject to spasticity may benefit from a larger Range setting, since it will increase the likelihood of a user successfully operating the switch (usually by a pass-through movement). It may help to introduce some Ignore Repeats time in this scenario, since this category of users will sometimes be unable to control the return of their hand to its resting position without making a second (return) pass through the beam.

Conversely, users with conditions such as muscular dystrophy or a motor neurone disease may benefit from a smaller range, allowing tiny movements to be recognised, such as a finger entering the sensor zone by just a few millimetres (a wave movement), or a cheek being puffed out (a top-down movement).

The information provided above represents just a small number of likely use-cases by way of example. Experimentation will lead to the correct position and method(s) of operation for each particular user.




Figure 2a – Pass-through activation




Figure 2b – Wave activation




Figure 2c – Top-down activation




Figure 2d – Activation and Release points

Effect of bright daylight/sunlight

Although the unit uses a very accurately controlled light source, it is not completely impervious to the effects of strong daylight and direct sunlight. The effect is to reduce the achievable range of the unit, and it is therefore more noticeable at larger Range settings. When exposed to very strong light you may find that the point at which the object is first recognised is less than the Range setting implies. For example, if the unit is set to Range 9 (500mm), you could find that the actual operating maximum is reduced to half that when operating in direct sunlight. The effect is much less noticeable at smaller Range settings, so it's best to use large Range settings only when necessary, particularly if the unit is likely to be subject to strong light.

Warranty

ToFFEE BT is warranted for two years against defects in manufacture or component failure. The unit is designed for use in domestic, educational, health and commercial applications. Use outside these areas is not authorised and will invalidate the warranty. Unauthorised repair or modification, mechanical abuse, immersion in any liquid or connection to incompatible equipment will also invalidate the warranty.

ToFFEE BT is not authorised for use with safety critical, mission critical or life-sustaining apparatus of any kind.

Maintenance

Your **ToFFEE BT** has no user serviceable parts. If repair becomes necessary, the unit should be returned to Pretorian Technologies or an authorised supplier.

ToFFEE BT contains a rechargeable lithium-ion battery that isn't user replaceable. Although the unit uses the very latest battery and charging technology, it may eventually need to be replaced, in which case it should be returned to Pretorian Technologies or an authorised supplier.

Disposal of batteries is often subject to local laws. Please consult your local authorities for information relating to your locale. Never dispose of a battery in fire.

Technical Specification

Ingress rating: IP30

Bluetooth Range: Approximately 15m (45 feet)

Mounting holes: 3mm threaded, 31.75mm pitch, 15mm maximum length

Troubleshooting Guide

Symptom	Possible Cause	Remedy
ToFFEE BT is not 'discoverable' on my device.	<ul style="list-style-type: none">• Battery is flat• Unit is in Sleep mode• Bluetooth is not enabled in device.• Unit is already connected to another device in range.	<ul style="list-style-type: none">• Charge the unit• Press Power button [7]• Enable Bluetooth.• Use 'forget this device' in Bluetooth menu of other device to make unit available again.
No power to ToFFEE BT	<ul style="list-style-type: none">• Battery is flat• Unit is in Sleep mode	<ul style="list-style-type: none">• Charge the unit• Wake unit up by pressing Power button [7].
ToFFEE BT does not wake up when I wave my hand in the beam	<ul style="list-style-type: none">• Unit has been put to sleep manually.	<ul style="list-style-type: none">• Wake unit up by pressing Power button [7].
Switch activates multiple times	<ul style="list-style-type: none">• Operating on edge of sensor zone• Operating at the farthest point of the sensor zone	<ul style="list-style-type: none">• Move hand further into the centre of the beam• Move hand closer or increase the Range setting <p>If the above measures don't help, consider using One-Shot mode, including some Ignore Repeats time, or both.</p>
ToFFEE BT operates randomly with no target in the sensitive zone	<ul style="list-style-type: none">• Cover glass [1] is dirty	<ul style="list-style-type: none">• Clean cover glass [1] with a soft, slightly damp cloth and polish.
When I select timed play, the music stops.	<ul style="list-style-type: none">• Music/video was already playing	<ul style="list-style-type: none">• Ensure playback is paused before selecting a timed play function.

Trademarks used in this document are for illustrative purposes only, are the property of their respective owners, and are acknowledged.



Product
Page



Unit 37 Corringham Road Industrial Estate
Gainsborough, Lincolnshire DN21 1QB UK
Tel +44 (0) 1427 678990 **Fax** +44 (0) 1427 678992



European Healthcare & Device Solutions Ltd.
Stratton House, Bishopstown Road, Cork T12 Y9TC, Republic of Ireland
(Do not use this address for correspondence with Pretorian).



S040226

Firmware revision is 164.0 onwards.

www.pretorianuk.com